

EASY-TO-USE FACE SHIELDS

- 1) .015 Clear PVC shield
- 2) Anti fog on interior side
- 3) Soft Form Forehead Open Cell Foam
- 4) 12" wide and 12 3/4 tall
- 5) Connects with PVC re-enforced Velcro Hand Straps
- 6) Individually bagged, boxes of 50

*Please read attached disclaimer of product use.

FROM OUR CUSTOMERS:

"I like these, they cover so much more, if you have to tilt your head upwards you are still covered..."

"Easy to clean and reuse."



FOR OTHER COVID-19 RESPONSE ITEMS LIKE LEVEL-1 GOWNS, CUSTOM PARTITIONS, REMOTE OFFICE BACKDROPS, CUSTOM BRANDED FABRIC MASKS, PLEASE CONTACT YOUR SYNERGY REPRESENTATIVE.

Alternate Disclaimer for PPE Products

Nomadic Display products have been manufactured and distributed during the COVID-19 Pandemic, when FDA cleared face masks are not readily available. Our products are not designed for use in medical settings when protective devices that have received FDA clearance are available and should not be considered an appropriate substitute for FDA cleared medical devices including surgical masks, when such devices are available for use. Our products are not recommended for use in a surgical setting or where the wearer may be exposed to hazardous bodily fluids. No claims are made as to our products' ability to prevent or reduce exposure, whether through dermal contact or respiration, to microbes, viruses, pathogens, or other agents that may cause harm to the wearer following contact or exposure. Our products are sold or made available for use only after the potential wearer has read the above statement. By using our products, the wearer acknowledges that they have received the above notice.

Nomadic Display makes no claim or warranties, real or implied, as to the effectiveness of our products in the prevention of infection of any type. Everyone is urged to follow CDC guidelines during this crisis. Please consult www.cdc.org for information and updates.

Labeling requirements for Level 1 & 2 Gowns and Masks

As of March 31st, 2020

According to the FDA, inside the box and on any paperwork is acceptable. A label does not need to be sewn in. The outer case in which it is shipped and distributed should have this info as well.

Contact Suzanne B. Schwartz, MD, MBA with further questions, not answered on the FDA website:
<https://www.fda.gov/medical-devices/general-hospital-devices-and-supplies/personal-protective-equipment-infection-control>

Suzanne B Schwartz, MD, MBA

Deputy Director (& Acting Office Direct) Office of Strategic Partnerships & Technology Innovation

Center for Devices and Radiological Health (CDRH)

Office of Strategic Partnerships and Technology Innovation (OST)

U.S. Food and Drug Administration

Wo66, Room 5410

Tel: 301-796-6937

Cell: 202-841-9996

Suzanne.Schwartz@fda.hhs.gov